

# AP-IgG Conjugation Kit

## For labeling 1 x 1 mg IgG

### **REAGENT STORAGE**

The kit is shipped on blue ice. Please store kit components as described below.

Kit Component	Storage Temp	Storage Notes
Concentrated Activator	-20°C	Keep the vial in the desiccated container as supplied in the kit
AP-Z™	-20°C or 2-8°C	Does not need to be kept desiccated.
Quenching Reagent	-20°C or 2-8°C	Does not need to be kept desiccated.

### **INTRODUCTION**

Alkaline Phosphatase is widely used as a fluorescent label in immunochemistry techniques such as flow cytometry and cellular analysis. Preparation of bright, stable, and reproducible antibody-AP conjugates is one of the biggest challenges for development of high-quality fluorescent reagents for cellular analysis and flow cytometry. The AP-IgG conjugation kit utilizes a novel chemistry to generate bright and highly reproducible AP-IgG conjugates with a simple procedure

### **FEATURES**

- Liquid-based reagents.
- Completely scalable: conjugate anywhere from 10 µg to 1-gram IgG per reaction.
- Supplies sufficient activated AP to conjugate all IgG at a 2:1 AP:IgG ratio.
- Highly efficient AP incorporation - purification not usually necessary.
- Customize the AP:IgG ratio to create optimized conjugates for different applications.
- Conjugates have greatly improved stability vs Lightning-Link™ and traditional chemistry.

### **PRODUCT and CONTENTS**

<b>Catalog Number</b>	<b>CNJ93103-001</b>
<b>For Labeling:</b>	<b>1 x 1 mg IgG</b>
Concentrated Activator	10 µL
AP-Z™ - Activated AP (7.4 mg/ml)	200 µL
1X Quenching Reagent	60 µL

### **ADDITIONAL REAGENTS REQUIRED BUT NOT SUPPLIED**

1X Phosphate Buffered Saline (1X PBS), pH 7.2-7.5

Deionized water (dH<sub>2</sub>O)

1.5 ml microcentrifuge tubes

## **SHELF LIFE**

The performance of the product is guaranteed for a minimum of 12 months when stored as directed.


## **IgG Requirements**

The IgG to be labeled should be at a minimum concentration of 0.8 mg/ml in pure 1X PBS and should not contain any preservatives or carriers such as sodium azide, Proclin 300 or BSA.

## **AP:IgG Molar Ratio**

This kit utilizes a 2:1 AP:IgG molar ratio which is optimal for most conjugations reaction. However, lower or higher ratios may give better results depending upon the antibody characteristics and the intended end-use. To change the AP:IgG molar ratio, vary the volume of AP-Z™ added to the conjugation reaction (Step 8).

## **CONJUGATION PROCEDURE - OVERVIEW**



## **BEFORE BEGINNING THE PROCEDURE**

Remove the Concentrated Activator from the freezer. **Important:** Allow sufficient time to let the container and contents come to room temperature before opening the outer and inner vials.

**Note:** The jar containing the Activator can be removed from the freezer up to 24 hours before use.

## **DETAILED CONJUGATION PROCEDURE**

1. Measure the absorbance of the IgG solution at 280 nm using PBS as a blank. Divide the A280 by 1.40 to obtain the IgG concentration in mg/ml.
2. Dilute IgG to 1.20 mg/ml in 1X PBS (0.80 – 1.4 mg/ml is acceptable).
3. Add 1 mg of IgG solution to a new microcentrifuge tube.
4. Prepare a working dilution (1X) of Activator from Concentrated Activator in deionized water:
  - a. Add 2.0  $\mu$ L of Concentrated Activator to 60  $\mu$ L of deionized water.
  - b. Immediately vortex to mix the solution thoroughly.

**Note:** The 1X Activator must be used within 5 minutes of preparation. If more than 5 minutes passes before use, discard the 1X Activator and prepare a fresh solution.

5. Add 10.0  $\mu$ L of 1X Activator to the 1 mg aliquot of IgG and then mix thoroughly by gentle vortexing.
6. Incubate the solution at room temperature for 1 hour.

**Note:** A longer incubation is not harmful and even overnight incubations will be successful.

7. Desalt the complete reaction volume into pure 1X PBS. We recommend using Pierce Zeba desalting spin columns with a 7 Kd MW cutoff for small volumes of IgG. Use of gravity desalting columns, dialysis, and extensive washing with centrifugal filter units for desalting is also acceptable.

**Note:** The activated IgG is stable and can be stored at 2-8°C for at least 4 months.

8. Add 200  $\mu$ L of AP-Z™ to the desalted, activated IgG and mix by gentle vortexing.
9. Incubate the solution at room temperature for 2-24 hours.

**Note:** Usable conjugates are produced after only 2 hours incubation. Larger and more potent conjugates will be produced after longer incubations.

10. Add 50  $\mu$ L of Quenching Reagent to the reaction and mix by gentle vortexing.
11. Incubate the solution at room temperature for 1 hour.

**Note:** A longer incubation is not harmful and overnight incubations will be successful.

12. Conjugate is ready for use. Store at 2-8°C.

**Note:** To improve conjugate performance, it may help to purify the conjugate from the unincorporated AP and reaction components by size exclusion chromatography.

## RECOMMENDED ACCESSORIES

For desalting IgG after activation - Order from ThermoFisher:

<b>Sample Size</b>	<b>Description</b>	<b>Cat #</b>
2 – 12 µL	Zeba Spin Desalting Columns, Micro (75µL), 7K MWCO	89877, 89878
30 - 130 µL	Zeba Spin Desalting Columns, 0.5mL, 7K MWCO	89882, 89883
200 – 700 µL	Zeba Spin Desalting Columns, 2mL, 7K MWCO	89889, 89890
500 – 2000 µL	Zeba Spin Desalting Columns, 5mL, 7K MWCO	89891, 89892
700 – 4000 µL	Zeba Spin Desalting Columns, 10mL, 7K MWCO	89893, 89894

For concentrating IgG before or after activation or for concentrating the final conjugate – Order from MilliporeSigma:

<b>Sample Size</b>	<b>Description</b>	<b>Cat #</b>
Up to 500 µL	Amicon Ultra-0.5 Centrifugal Filter Unit with Ultracel-50 membrane	Z740176
Up to 2 mL	Amicon Ultra-2 Centrifugal Filter Unit with Ultracel-50 membrane	UFC205024
Up to 4 mL	Amicon Ultra-4 Centrifugal Filter Unit with Ultracel-50 membrane	UFC805008
Up to 15 mL	Amicon Ultra-15 Centrifugal Filter Unit with Ultracel-50 membrane	Z648000