

Langford Locks Kidlington Oxfordshire OX5 1LH

Telephone: +44(0)1865 595230

Antibody Datasheet

Name: Goat Anti-Salmonella CSA-1 Antibody

Product Code: 5310-0322

Batch #:

Date of Manufacture:

Product Description: Affinity purified polyclonal antibody to the Salmonella Common Structural Antigens (CSA) made in goat (BacTrace®). Isolated from a serum pool of goats immunized with different serotypes of Salmonella. This antibody is broadly reactive to Salmonella and recognizes most Salmonella serotypes.

Amount: 1.0 mg (as determined by UV absorbance at 280 nm)

Presentation: Lyophilized

Rehydration: Two methods for rehydration and storage are recommended to meet most needs, Procedure A, using 50% glycerol, is preferred. It eliminates freezing at -20 °C and is an effective biological inhibitor. At a working dilution, the level of glycerol is too small to affect most assays.

Procedure A: 50% Glycerol

- o Solution Preparation: Mix 1 part glycerol with 1 part reagent quality water.
- Rehydration: Transfer 1 mL of the 50% glycerol solution to the product vial. Rotate the vial until the lyophilized pellet is totally dissolved. Antibody concentration will then be 1.0 mg/mL.
- Storage: The rehydrated product is stable for a minimum of 1 year at 2-8 °C.

Procedure B: Acetic Acid/Bi-Carbonate

- Solution Preparation: Prepare immediately before use:
 - o M Acetic Acid Solution: Mix 25 μL of Glacial Acetic Acid and 40 mL reagent quality water.
- 0.177M Carbonate-Bicarbonate Solution: Dissolve 1.09 g of Na2CO3 and 0.63 g of NaHCO3, in
 100 mL reagent quality water.
- o Rehydration: Add 0.1 mL of 0.01M Acetic Acid to 1 mg of antibody. Rotate vial until totally dissolved. Add 0.1 mL of the Carbonate-Bicarbonate Solution and rapidly mix to obtain a clear or opalescent solution. Heat at 37°C for 30 minutes and allow to slowly cool to room temperature. Final concentration is 5.0 mg/mL.
- Storage: This product may be stored for up to one week at 2-8°C; thereafter, it should be stored at -20 °C. At -20 °C, product is stable fur a minimum of 1 year.





www.thenativeantigencompany.com VAT No: 102038475

NAC.Contact@LGCGroup.com Registration No. 7386339



Langford Locks Kidlington Oxfordshire OX5 1LH

Telephone: +44(0)1865 595230

Stabilizer: No stabilizers added.

Preservative: No preservatives added. Non-sterile.

Purification: Affinity purified

Specificity: This product is designed to react broadly with Salmonella enterica. This includes serotypes of group A, B, C, D, E, and other antigen classes. The intensity of reaction may vary depending on the strain of Salmonella tested. Product is tested via ELISA techniques as applicable.

Cross-Reactivity: Antibody may also show some crossreactivity to related Enterobacteriaceae.

Application Notes: Optimal working concentrations should be determined experimentally. Prepare working dilution in PBS or other buffer. Dilution is not recommended for long-term storage. Suggested dilutions are listed below. In many cases, the antibody may be diluted further than indicated:

Capture Antibody: 10 µg/mL Primary Antibody: 1-2 µg/mL

Usage Guidelines

Storage: +4°C

Storage Guidelines: Store $+4^{\circ}$ C until rehydrated. Stable for a minimum of 1 year when stored at $+4^{\circ}$ C. Avoid multiple freeze-thaw cycles.

X			
QC			
Χ			
QA			

Products are for Research Use or for Further Manufacturing Use only. Not for Diagnostic or Therapeutic Use.



